

BARIATRIC WHEELCHAIR MOVER MBWC SPECIFICATIONS



Australia Wide Sales & Service

sales@materialshandling.com.au | www.materialshandling.com.au | 1300 65 00 35

MATERIALS
Handling PTY LTD *Working with ease.*

BARIATRIC WHEELCHAIR MOVER MBWC SPECIFICATIONS

WHEELCHAIR Mover^{AP}



Can your staff move a 250kg passenger and wheelchair across a carpet or up jetways smoothly, safely and effortlessly?

Introducing the WheelChair Mover– the industry’s first power-assist device making passenger wheelchair transport safe, smooth and easy. Designed by safety experts for use by airport personnel, the WheelChair Mover fits all standard wheelchair sizes and models in your terminal, easily transporting passengers safely and securely thru hallways, in elevators, and up jetways while reducing potential injury or physical strain to personnel. The WheelChair Mover’s quiet, smooth movement means a gentle ride while it’s maneuverability and ease of use means personnel move safely and efficiently.

Features

- Variable speed for passenger transport comfort
- Positive wheelchair locking mechanism
- Compact design for maneuverability
- Easy to use controls with system status indicators
- Floor friendly non-marring solid rubber tires
- Automatic brake prevents rollback
- Charger interlock prevents drive away while batteries are on charge
- Charging cord with convenient storage area

Specifications

- Dimensions: 610 W x 712 L x 1080 H mm 150 kg
- Electrical: 240 VAC 50 Hz 0.75 A
- Motor: 0.5 hp, 24 V, geared front wheel drive
- Batteries: 2 x 12 VDC sealed lead acid maintenance free
- 6 Amp multi-stage Smart Charger with industrial cord
- Capacity: 250 kg patient plus wheelchair
- Adapts to any standard wheelchairs: 450-610 mm wide



Powerful and safe



Easy to use controls



Manoeuvrable

Australia Wide Sales & Service

sales@materialshandling.com.au | www.materialshandling.com.au | 1300 65 00 35

MATERIALS
Handling PTY LTD Working with ease.